HemaSphere

Comment Open Access



TP53 Mutations and Clinical Outcome in Chronic Lymphocytic Leukemia: Is a Threshold Still Needed?

Riccardo Bomben¹, Antonella Zucchetto¹, Federico Pozzo¹, Erika Tissino¹, Tamara Bittolo¹, Jacopo Olivieri², Annalisa Chiarenza³, Francesco Zaja⁴, Maria Ilaria Del Principe⁵, Davide Rossi^{6,7}, Valter Gattei¹

Correspondence: Valter Gattei (vgattei@cro.it).

n a recent issue of *HemaSphere*, Pandzic and colleagues¹ provide evidence that the use of next generation sequencing (NGS) in chronic lymphocytic leukemia (CLL) clinical practice is equally robust in detecting variants with allele frequency (VAF) between 5% and 10% respect to those with VAF \geq 10%, and suggest to harmonize NGS methodologies across laboratories and to lower the canonical cutoff from \geq 10% VAF.

This is an important issue to address, since it is now emerging from the scientific literature that the detection of TP53 mutations also at the subclonal level, that is, far below the 10% VAF cutoff, may have a clinical impact and should be therefore annotated for a more correct clinical management of CLL patients. In fact, since the first seminal observation,3 other more recent retrospective studies based on ultra-deep NGS have confirmed that TP53 mutations can be present in tumor cell populations at very low-VAF,⁴⁻⁷ and that the presence of these extremely small TP53 clones, even with VAF <1.0%, was associated with inferior survival in the context of chemoimmune therapy.5-7 Recently, it was reported that also in the setting of therapies with the novel BCR inhibitor ibrutinib, TP53 mutations/disruptions were independent negative prognostic factors,^{4,8} although the impact of low-VAF TP53 mutations upon treatment with targeted agents remains currently unexplored.

In general, all the above-mentioned studies,³⁻⁷ rather than identifying a threshold of clinical relevance, seem to point toward a "yes/no effect" of *TP53* mutations on CLL outcome. Biologically speaking, this is in line with the basic principle of fitness where the fittest clones produce viable offspring capable

HemaSphere (2023) 7:4(e855).

http://dx.doi.org/10.1097/HS9.00000000000855.

of surviving and reproducing and, in the presence of selective pressures from the environment, for example, due to chemotherapy, eventually taking over the *TP53* wild-type cell population.³

Nevertheless, the European Research Initiative on CLL (ERIC) recommend that only pathogenic/likely pathogenic *TP53* variants with VAF >10% should be reported,² a cautionary approach proposed to avoid not only the risk of possible technical artifacts in the identification of very low-VAF *TP53* mutations, but also for the lack of solid evidence and/or prospective studies on the clinical impact of low-VAF TP53 variants. Recently, however, ERIC has launched an interlaboratory comparison of NGS methods to detect low-VAF variants; the multilaboratory initiative aimed to harmonize methods detecting TP53 variants below 10% VAF to be applied for the multicentric study analyzing their prognostic/predictive impact in CLL.⁹

Pandzic and colleagues¹ thoroughly and brilliantly discussed all the limitations, drawbacks, and pitfalls of the use of the 10% VAF threshold, by also highlighting some paradoxical situations that may happen in the clinical practice for a different evaluation of the *TP53* mutational burden due to the presence or not of a parallel deletion of the contralateral 17p allele. Yet, they provide extra-assay evidence (ie, by droplet digital polymerase chain reaction [PCR]) against a lowering of the threshold below the 5% VAF.¹

Here we summarize the internal recommendations we are used to follow for NGS evaluation of *TP53* mutational status in CLL, including wet laboratory and bioinformatics, as reported in a recent publication of ours,⁷ where we described the set-up of an ultra-deep NGS procedure capable to successfully detect mutations below the 1% VAF. A custom script of the bioinformatics pipeline, as described herein and in reference 7, can be retrieved at website address https://github.com/gamabunta313/ TP_SNP-calling/. More details are reported in reference 7.

Briefly, analysis of *TP53* mutations is performed with an amplicon-based strategy, covering exons 2-11, in keeping with ERIC recommendations.² Amplicon libraries are generated starting from 40 ng of DNA (~6000 diploid genomes), a quantity capable to successfully detect mutations below the 1% VAF in the context of our procedures.⁷ Multiplex PCR products are generated using a high-fidelity DNA polymerase and tagged with specific indices and paired-end sequenced in an Illumina MiSeq instrument.

FASTQ files are aligned to the hg19 reference with Burrows-Wheeler Aligner-MEM algorithm, and allele variants are called by FreeBayes with nonstringent parameters. The overall

¹Clinical and Experimental Onco-Hematology Unit, Centro di Riferimento

Oncologico di Aviano (CRO), IRCCS, Aviano (PN), Italy

²Clinica Ematologica, Centro Trapianti e Terapie Cellulari "Carlo Melzi" DISM, Azienda Ospedaliera Universitaria S. Maria Misericordia, Udine, Italy

Azienua Ospeualiera Oniversitaria S. Mana Misencorula, Ouline, ita

³Division of Haematology, Ferrarotto Hospital, Catania, Italy

⁴Department of Medical, Surgical and Health Sciences, University of Trieste, Italy ⁵Division of Haematology, University of Tor Vergata, Rome, Italy ⁶Haematology, Institute of Oncology Research, Bellinzona, Switzerland ⁷Oncology Institute of Southern Switzerland, Bellinzona, Switzerland ⁸Oncology Institute of Southern Switzerland ⁸Oncology Institute of Southern Switzerland, Bellinzona, Switzerland ⁹Oncology Institute of Southern Switzerland ⁹Oncology Institute of Southe




Figure 1. Overall survival of patients stratified according to the *TP53* mutation load. CLL cases from a retrospective cohort originally published in⁷ were reclassified according to the standard cutoff for *TP53* mutation evaluation ($\geq 10\%$ VAF), to the cutoff suggested by¹ ($\geq 5\%$ VAF), and to the <5\% VAF cutoff, as discussed in the present commentary. (A) Kaplan-Meier curves comparing OS computed from the date of *TP53* mutation evaluation⁷ of CLL 1220 patients stratified into 4 classes according to the *TP53* mutational load: *TP53* wild-type cases (green line, n = 1052), cases with <5\% VAF of *TP53* mutation (blue line, n = 64), cases with VAF of *TP53* mutations between 5.0% and 10.0% (purple line, n = 12), and cases with $\geq 10.0\%$ VAF of *TP53* mutation (red line, n = 92). (B) Kaplan-Meier curves comparing OS computed from the date of treated CLL patients stratified into 4 classes according to the *TP53* mutation to due of treatment⁷ of 552 treated CLL patients stratified into 4 classes according to the *TP53* mutational load: *TP53* wild-type cases (green line, n = 49), cases with <5% VAF of *TP53* mutation (blue line, n = 38), cases with VAF of *TP53* mutations between 5.0% and 10.0% (VAF of *TP53* mutation (blue line, n = 38), cases with VAF of *TP53* mutations between 5.0% and 10.0% (purple line, n = 4), and cases with $\geq 10.0\%$ VAF of *TP53* mutation (blue line, n = 61). *P* values refer to the log-rank test. CLL = chronic lymphocytic leukemia; 0S = overall survival; VAF = variants with allele frequency.

strategy is designed to yield a minimum coverage of 2000X in 100% of the analyzed positions for each sequence, a prerequisite for $\leq 1\%$ variant calling. To calculate random/systematic errors, a database of *TP53* wild-type cases is utilized (362 cases in reference 7). A *TP53* mutation is accepted as valid if (the following conditions have to be both fulfilled): (i) presented a VAF outdistancing for at least 2.75 standard deviations the mean of the transformed VAF at the specific *TP53* position as retrieved in the wt cases; (ii) validated by Fisher exact test after Bonferroni correction (P < 0.01). The minimal allelic fraction found in reference 7 was 0.3%.

Even if robust bioinformatic workflows are employed to call true variants from background error noise, to exclude possible random polymerase errors despite the use of a high fidelity enzyme, it would be advisable to reanalyze mutation less than 2.0% of VAF. In addition, the detected variants should be checked using locus-specific database, that is, the TP53 Database (https://tp53.isb-cgc.org/), for interpretation of *TP53* variants to discriminate between possible neutral and/or pathogenetic variants; programs like Mutalyzer (https://mutalyzer. nl/) and Sequence Variant Nomenclature (http://varnomen.hgvs. org/) are used to report the variant nomenclature according to the guidelines of the Human Genome Variation Society.

By such an approach, and re-classifying the retrospective cohort investigated in reference 7 (n = 1220) according to the cutoff suggested by 1, no difference was found in term of overall survival (OS) when comparing CLL cases with VAF lower than 5%, between 5% and 10% and upper than 10%, all the 3 categories experiencing significantly shorter OS intervals than *TP53* wild-type CLL (Figure 1A). The same held true when the analysis was circumscribed to the 552 treated CLL cases, in which OS was computed from treatment initiation (Figure 1B).

In conclusion, it is time to go beyond the 10% VAF cutoff for *TP53* mutation evaluation in CLL by lowering the cutoff at least at the value of 5% VAF, as suggested by reference 1, or even below, when the technical quality is guaranteed by a solid laboratory and bioinformatics procedure; in this regard, some suggestions from reference 7 are summarized in the present commentary. However, it is to remind that any novel cutoff for *TP53* mutation evaluation should be upfront validated by ad-hoc trials and cannot be introduced in the clinical practice in the absence of such a validation. Finally, *TP53* variants analysis through a very sensitive NGS approaches should be centralized in few recognized/certified laboratories to fulfill all the criteria reported above and to guarantee the same standard of precision medicine for each patients.

AUTHOR CONTRIBUTIONS

All authors conceptualized, wrote, and finalized the manuscript.

DISCLOSURES

The authors have no conflicts of interest to disclose.

SOURCES OF FUNDING

The authors declare no sources of funding.

REFERENCES

- Pandzic T, Ladenvall C, Engvall M, et al. Five percent variant allele frequency is a reliable reporting threshold for TP53 variants detected by next generation sequencing in chronic lymphocytic leukemia in the clinical setting. *HemaSphere*. 2022;6:e761.
- Malcikova J, Tausch E, Rossi D, et al. ERIC recommendations for TP53 mutation analysis in chronic lymphocytic leukemia - update on methodological approaches and results interpretation. *Leukemia*. 2018;32:1070–1080.

HemaSphere (2023) 7:4

- Rossi D, Khiabanian H, Spina V, et al. Clinical impact of small TP53 mutated subclones in chronic lymphocytic leukemia. *Blood*. 2014;123:2139–2147.
- 4. Brieghel C, Aarup K, Torp MH, et al. Clinical outcomes in patients with multi-hit TP53 chronic lymphocytic leukemia treated with ibrutinib. *Clin Cancer Res.* 2021;27:4531–4538.
- Malcikova J, Pavlova S, Kunt Vonkova B, et al. Low-burden TP53 mutations in CLL: clinical impact and clonal evolution within the context of different treatment options. *Blood*. 2021;138:2670–2685.
- Nadeu F, Delgado J, Royo C, et al. Clinical impact of clonal and subclonal TP53, SF3B1, BIRC3, NOTCH1, and ATM mutations in chronic lymphocytic leukemia. *Blood.* 2016;127:2122–2130.
- Bomben R, Rossi FM, Vit F, et al. TP53 mutations with low variant allele frequency predict short survival in Chronic Lymphocytic Leukemia. *Clin Cancer Res.* 2021;27:5566–5575.
- Morabito F, Tripepi G, Del Poeta G, et al. Assessment of the 4-factor score: Retrospective analysis of 586 CLL patients receiving ibrutinib. A campus CLL study [Published online 2021]. Am J Hematol. 2021;96:E168–E171.
- 9. Pavlova S, Malcikova J, Radova L, et al. Interlaboratory comparison of NGS methods for detection of TP53 Variants <10% VAF: the first phase of an ERIC multicenter study on the clinical impact of low-frequent TP53 variants in CLL. *HemaSphere*. 2020;4(S1):285–286.